

SpiSOP update 2.3.8.3

New standalone **SpiSOP version 2.3.8.3** is now online

- Improved browser for sleep scoring:
 - Quickstart to load EDF or Brainvision files immediately without setup of SpiSOP
 - Automated sleep scoring service (beta support) from Z3Score (based on EEG and EOG activity), also see their introduction video or contact them on testdrive@z3score.com
 - ...and “skip high confidence epochs” while (re-)scoring the hypnogram.
 - Comparison of up to 4 hypnograms that were (automatically) scored
 - Supporting >100 channels to be viewed better
 - etc.
- Added new Rapid Eye-movement (REM) detection (function name `remsMaAd`) with great help from Marek Adamczyk in using his code and methodology. be aware, it is not extensively tested in SpiSOP yet.
- Output is now created automatically in respectively named subfolders ordered per analysis/run name and function name.
- Many other small feature additions and bug fixes, e.g fixing of reading EDF that have inversed physical maximum and minimum (formerly this created errors or wrong values/analyses).